

RECORDS ADMINISTRATION  
R0080950

DP-MS-67-114

671491

IMPROVED ACTIVATION ANALYSIS OF TRACE  
ELEMENTS IN FRESH WATER SOLIDS\*

by

R. F. Overman

Savannah River Laboratory  
E. I. du Pont de Nemours and Co.  
Aiken, South Carolina 29801

SRL  
RECORD COPY

Proposed for publication in Limnology and Oceanography.

\*The information contained in this article was developed during the course of work under Contract AT(07-2)-1 with the U. S. Atomic Energy Commission.

5/13/68

This document was prepared in conjunction with work accomplished under Contract No. AT(07-2)-1 with the U.S. Department of Energy.

#### DISCLAIMER

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

This report has been reproduced directly from the best available copy.

Available for sale to the public, in paper, from: U.S. Department of Commerce, National Technical Information Service, 5285 Port Royal Road, Springfield, VA 22161

phone: (800) 553-6847

fax: (703) 605-6900

email: [orders@ntis.fedworld.gov](mailto:orders@ntis.fedworld.gov)

online ordering: <http://www.ntis.gov/help/index.asp>

Available electronically at <http://www.osti.gov/bridge>

Available for a processing fee to U.S. Department of Energy and its contractors, in paper, from: U.S. Department of Energy, Office of Scientific and Technical Information, P.O. Box 62, Oak Ridge, TN 37831-0062

phone: (865)576-8401

fax: (865)576-5728

email: [reports@adonis.osti.gov](mailto:reports@adonis.osti.gov)

IMPROVED ACTIVATION ANALYSIS OF TRACE  
ELEMENTS IN FRESH WATER SOLIDS

## ABSTRACT

Trace elements in water or solids can be easily determined after neutron activation with the high-resolution lithium-drifted germanium detector and gamma-gamma coincidence equipment. The gamma energies are determined more precisely by the germanium detector than is possible with sodium iodide detectors. Observed activities are identified by comparing them with the known decay characteristics of radioactive nuclides. The calculations do not require computers, and the number of chemical separations is reduced. This technique was used to determine, without chemical separation, sodium, bromine, arsenic, chromium, iron, cobalt, scandium, lanthanum, and zinc in particulate and evaporated filtrate limnological samples.

## INTRODUCTION

Trace elements can be accurately determined with neutron activation analysis. Methods for analyzing marine organisms have been reviewed by Fukai and Meinke (1959 and 1962). Methods for analyzing water for trace metals have been compared by Hume (1967). Neutron activation services, for a simple irradiation to the complete determination of complex samples, are now available commercially; 21 companies offering such services have been listed by Lyon and Miller (1967).

Most activation methods quantitatively determine the gammas emitted by the nuclides of interest. Interferences in complex gamma spectra can be minimized by electronically or mathematically subtracting known spectra from the sample spectrum -- denoted "spectrum stripping" -- and by chemically separating the mixture of nuclides. Both methods assume that the components of the mixture are known; neither is ideal for all mixtures. Activation techniques may become more applicable either by improving the chemical separation techniques or by improving the resolution of the gamma detection system. The high resolution of the lithium-drifted germanium detector improves the resolution.

This paper demonstrates the advantages of lithium-drifted germanium, Ge(Li), detectors, both alone and in conjunction with sodium iodide, NaI(Tl), detectors for gamma-gamma coincidence measurements. A description of the Ge(Li) detector system is given by Camp (1967). The high resolution (4 kev or better) of these detectors permits resolution of most of the gammas in activated samples; the energy and count rate of each photopeak can

be easily calculated. With gamma-gamma coincidence techniques, the observed gamma energies can be assigned to particular nuclides, and absolute counting is possible (Siegbahn 1965).

This paper presents techniques for determining trace elements in solid samples obtained from a fresh water pond. Sodium, bromine, arsenic, iron, cobalt, chromium, scandium, lanthanum, and zinc were determined simultaneously without chemical separation.

## EXPERIMENTAL

### Instruments and Detectors

Ge(Li) and NaI(Tl) gamma detectors were used. The Ge(Li) detector was manufactured by Isotopes, Inc. It is an 8 cm<sup>2</sup> by 10-mm-thick planar detector with a 3.7 kev resolution for the 1.3 Mev peak. The detector is cooled by liquid nitrogen; however, the first-stage field effect transistors in the preamplifier are at room temperature. The resolution could be increased by cooling the transistors, as reported by Cline (1967). The two detectors used for the gamma-gamma coincidence measurements were 3-inch-diameter by 3-inch-long NaI(Tl).

A 4096 channel analyzer, Nuclear Data, Inc., with a dual 4096 channel analog-to-digital converter (ADC) was used. A Tennelec TC-200 amplifier was used with the Ge(Li) detector. The dual ADC presented the coincidence data from both NaI(Tl) detectors to the analyzer; thus all coincidence data were available at one time.

### Sample Preparation and Irradiation

The samples were solids filtered or evaporated from a fresh water pond. They were ashed to remove organic matter, weighed,

and wrapped in aluminum foil for irradiation. Pieces of the aluminum foil were also irradiated as blanks because the solids were not removed from the foil before counting. Some bromine or arsenic activity may volatilize in ashing; however, this possibility was not investigated.

The samples were placed in a water-tight aluminum capsule and irradiated in a reactor to about  $10^{17}$  neutrons/cm<sup>2</sup>. Exposures  $>10^{17}$  neutrons/cm<sup>2</sup> would be required to increase the sensitivity of the method so that elements with half-lives greater than 48 hours could be determined. Because of the high <sup>24</sup>Na activity induced in the aluminum capsule, the samples decayed for six days before they were counted.

#### Observed Gamma Spectra

The energies and assigned nuclides of the observed gamma peaks are listed in Table 1. The nuclide assignments were based on published decay schemes (Lederer, Hollander, and Perlman 1967), half-lives, and known occurrence in nature. The peaks in low abundance were not assigned.

The gamma spectra of an evaporated filtrate, which had decayed seven days, are shown in Figure 1. The two spectra illustrate the advantage of the 3.7-keV resolution of Ge(Li) over the 60-keV resolution of NaI(Tl) detectors. Even after seven days of decay <sup>24</sup>Na and <sup>82</sup>Br are still predominant. The <sup>82</sup>Br peak at 0.554 MeV is about the same height as the <sup>82</sup>Br peak at 0.777 MeV. However, the observed peak at 0.554 MeV was higher than that at 0.777 MeV which indicated the presence of an additional nuclide. The unresolved nuclide was deduced to be <sup>76</sup>As.

A gamma-gamma coincidence spectrum of the evaporated filtrate sample used in Figure 1 was made six and seven days after reactor discharge. The spectrum in Figure 2 was taken six days after reactor discharge and shows that the 1.37 Mev gamma is in coincidence with a 2.75 Mev gamma. The decay scheme, shown in Figure 3, indicates that  $^{24}\text{Na}$  has these decay characteristics. The coincidence spectrum taken after seven days of decay, Figure 4 A and B, shows that the 0.777 Mev gamma is in coincidence with 0.554, 1.045, and 1.320 Mev gammas, and also the 0.554 Mev gamma is in coincidence with the 1.320 Mev gamma. The decay scheme in Figure 5 indicates that  $^{82}\text{Br}$  has these decay characteristics. Other nuclides are identified by gamma energies, half-lives, and the known natural occurrence of precursors in the activation reaction.

The gamma spectra of a particulate sample are shown in Figure 6. The coincidence spectrum after seven days of decay is similar to Figure 4 because  $^{82}\text{Br}$  is still a major contributor. The coincidence of 0.487 and 1.596 Mev gammas confirmed the presence of  $^{140}\text{La}$ . The coincidence spectrum after twenty days decay, Figure 7 A and B, shows that the 0.889 Mev gammas are in coincidence with the 1.120 Mev gammas, and the 1.173 Mev gammas are in coincidence with 1.332 Mev gammas. These coincidence gammas are assigned to  $^{46}\text{Sc}$  and  $^{60}\text{Co}$ , respectively.

## CALCULATIONS AND RESULTS

### Calculations from Counting Data

The disintegration rate of each nuclide was calculated from the Ge(Li) data. Each peak was taken to be 20 kev wide, and its baseline was taken to extend 10 kev on either side of the peak.

The count rates for the peak areas were corrected for the baseline count rate. The counting efficiencies for the various energies of gammas were determined by counting calibrated  $^{110m}\text{Ag}$  and  $^{60}\text{Co}$  sources. The branching ratios, half-lives, and abundances of the gammas for each nuclide were taken from Lederer, Hollander, and Perlman (1967). In the case of  $^{76}\text{As}$ , the  $^{82}\text{Br}$  contribution to the 0.554 Mev peak was calculated, and  $^{76}\text{As}$  was reported as the difference. The  $^{65}\text{Zn}$  concentration in Figure 6 was determined in a similar fashion.

#### Irradiation Calculations

The method of calculating the grams of naturally occurring nuclides in the sample is given in "Activation Analysis Handbook" by Koch (1960). The thermal cross sections were also taken from the Handbook. Since this irradiation was intended to determine the nuclides that were detectable, comparative standards were not run as they would be if the methods were put into routine use. The neutron exposure was determined by irradiating a known amount of cobalt along with the samples.

#### Results of Calculations

The results for ten samples, five particulate and five evaporated filtrate, are given in Table 2. (The samples are from different depths, and therefore not replicates.) The concentration is based on the weight of the solid sample after ignition and not on the amount of water filtered. The accuracy or precision of the method were not determined.

## SUMMARY AND CONCLUSIONS

This study shows ways in which the usual neutron activation analysis methods can be improved. Chemical separations for many elements can be omitted if a high resolution Ge(Li) detector and an analyzer with a large number of channels are used. Nuclides can be assigned to observed gamma peaks easier by using a multi-parameter analyzer and two gamma detectors to determine which gammas are in coincidence. These improvements should reduce significantly the number of man-hours required per analysis.

# REFERENCES

- Camp, D. C. 1967. Applications and Optimization of the Lithium - Drifted Germanium Detector System. Tech. Rept. UCRL 50156, Lawrence Radiation Laboratory, Livermore, Calif.
- Cline, J. E. and R. L. Heath. 1967. Annual Progress Report. Gamma Ray Spectrometry of Neutron-Deficient Isotopes. Tech. Rept. IDO-17222, Phillips Petroleum Co., Idaho Falls, Idaho.
- Fukai, Rinnosuke and W. W. Meinke. 1959. Trace analysis of marine organisms: A comparison of activation analysis and conventional methods. Limnol. Oceanog., 4: 398-408.
- \_\_\_\_\_. 1962. Activation analyses of vanadium, arsenic, molybdenum, tungsten, rhenium and gold in marine organisms. Limnol. Oceanog., 7: 186-200.
- Hume, D. N. 1967. Analysis of water for trace metals - present capabilities and limitations. In Advances in chemistry series, No. 67, Equilibrium concepts in natural water systems. American Chemical Society.
- Koch, R. C. 1960. Activation analysis handbook. Academic Press, New York.
- Lederer, C. M., J. M. Hollander, and I. Perlman. 1967. Table of isotopes. Sixth Edition, John Wiley and Sons, Inc., New York.
- Lyon, W. S. and F. J. Miller. 1967. Forensic applications of neutron activation analysis, isotopes and radiation technology. v. 4, No. 4, p. 325-340, USAEC Clearing House for Federal and Scientific Technical Information.
- Siegbahn, Kai [Ed.]. 1965. Alpha-, beta-, and gamma-ray spectroscopy. p. 539-555, North-Holland Publishing Company, Amsterdam.

TABLE 1. Observed Gamma Peaks and Assigned Nuclides

<u>kev</u>	<u>nuclide</u>	<u>kev</u>	<u>nuclide</u>
312		828	$^{82}\text{Br}$
320*	$^{51}\text{Cr}$	888*	$^{46}\text{Sc}$
330	$^{140}\text{La}$	1042	$^{82}\text{Br}$
336		1098*	$^{59}\text{Fe}$
347		1116*	$^{65}\text{Zn} + ^{46}\text{Sc}$
412		1171	$^{60}\text{Co}$
487*	$^{140}\text{La}$	1290	$^{59}\text{Fe}$
496		1317	$^{82}\text{Br}$
554*	$^{82}\text{Br} + ^{76}\text{As}$	1331*	$^{60}\text{Co}$
620	$^{82}\text{Br}$	1369*	$^{24}\text{Na}$
697	$^{82}\text{Br}$	1475	$^{82}\text{Br}$
776*	$^{82}\text{Br}$	1593	$^{140}\text{La}$
815	$^{82}\text{Br}$	1730	$^{24}\text{Na}$ (escape peak)

---

\* Peaks used in the calculations for the given nuclides.

TABLE 2. Concentration of Elements in Ashed Samples

Particulate	Co		Fe		Cr		La		Sc	
	c/m/mg	$\frac{\text{g/g sample}}{(10^{-6})}$	c/m/mg	$\frac{\text{g/g sample}}{(10^{-2})}$	c/m/mg	$\frac{\text{g/g sample}}{(10^{-4})}$	c/m/mg	$\frac{\text{g/g sample}}{(10^{-4})}$	c/m/mg	$\frac{\text{g/g sample}}{(10^{-4})}$
1	0.26	14	1.1	2.2	2.3	2.9	3.7	3.2	2.0	4.0
2	1.11	60	1.4	2.8	2.3	2.9	2.8	2.4	1.5	2.9
3	0.29	15	1.8	3.8	2.3	2.9	1.0	0.8	1.3	2.6
4	0.50	27	3.3	6.7	2.7	3.4	3.6	3.1	1.3	2.6
5	0.18	10	7.4	15.0	1.6	2.0	1.7	1.5	0.7	1.4
Evaporated Filtrate										
1	0.46	25	0.11	0.23	0.38	0.49	-	-	-	-
2	0.01	0.6	0.13	0.26	0.22	0.28	-	-	-	-
3	0.08	4	0.77	1.6	0.29	0.36	-	-	-	-
4	0.06	3	0.06	0.11	0.18	0.23	-	-	-	-
5	-	-	0.09	0.18	0.20	0.25	-	-	-	-

Particulate	Br		As		Na		Zn	
	c/mg	g/g sample (10 <sup>-5</sup> )	c/mg	g/g sample (10 <sup>-5</sup> )	c/mg	g/g sample (10 <sup>-2</sup> )	c/mg	g/g sample (10 <sup>-4</sup> )
1	1.7	3.8	0.50	0.68	6.6	0.67	0.24	4.0
2	2.1	4.7	0.23	0.31	8.8	0.89	0.22	3.7
3	0.6	1.4	3.8	5.1	12.0	1.2	0.20	3.3
4	3.2	7.4	1.0	1.4	10.6	1.1	0.08	1.3
5	2.4	5.5	2.7	3.6	3.0	0.3	0.38	6.3
Evaporated Filtrate								
1	22	52	12	16	56	5.7	0.06	1.0
2	23	53	10	13	56	5.7	0.01	0.2
3	12	28	5	7	25	2.5	0.05	0.8
4	21	47	9	12	49	4.9	0.03	0.5
5	18	40	19	26	42	4.3	-	-






FIG. 1 GAMMA SCANS OF EVAPORATED FILTRATE SEVEN DAYS AFTER DECAY COMPARING NaI(Tl) AND Ge(Li) DETECTORS



A



B

FIG. 2 GAMMA-GAMMA COINCIDENCE SPECTRUM OF EVAPORATED FILTRATE SIX DAYS AFTER DECAY




FIG. 3 DECAY SCHEME OF  $^{24}\text{Na}$



A



B

FIG. 4 GAMMA-GAMMA COINCIDENCE SPECTRUM OF EVAPORATED FILTRATE SEVEN DAYS AFTER DECAY




FIG. 5 DECAY SCHEME OF  $^{82}\text{Br}$






FIG. 6 GAMMA SCAN OF PARTICULATE SAMPLE SEVEN DAYS AFTER DECAY COMPARING NaI(Tl) AND Ge(Li) DETECTORS



A



B

FIG. 7 GAMMA-GAMMA COINCIDENCE SPECTRUM OF PARTICULATE SAMPLE TWENTY DAYS AFTER DECAY